


COLON ADENOCARCINOMA

RECOMMENDED MOLECULAR TESTING



	Test/CPT	Turnaround Time	Specimen Requirements
PCR	BRAF V600 81210 88381 may apply	Within 3-6 days of receipt (T&F)	• Formalin fixed, parafin-embedded (FFPE) tissue block or 10 FFPE tissue slides
	KRAS-IVD Exon 2 81275 88381 may apply	Within 3-6 days of receipt (T&F)	• Formalin-fixed, paraffin-embedded (FFPE) tissue block (preferred) or unstained slides (10 slides are usually sufficient)
	KRAS Exons 3,4 81276 88381 may apply	Within 5-8 days of receipt (T&Th)	
	NRAS Exons 2,3,4 81311 88381 may apply	Within 5-8 days of receipt (T&Th)	

*CPTs for reference only and may vary

Updated: August 2018